Description / composition

calc-i-oss™ is a bone graft substitute made from phase-pure beta-tricalcium phosphate (β -Ca₃(PO₄)₂). It is available as granules.

calc-i-oss™ is biocompatible and osteoconductive. In the defect, the granules form a macroporous structure that fosters formation of autologous tissue. The granule's high microporosity allows for taking up considerable amounts of blood. With the uptake of blood, factors that are beneficial for healing are made available.

calc-i-oss™ is fully resorbable.

calc-i-oss™ is synthetic and does not contain animal or human-derived substances.

calc-i-oss™ is sintered at over 1000 °C. The granulation and sinter process produces round, porous granules of high stability with only marginal inclination towards fine particle attrition.

Indications

calc-i-oss™ is a bone graft substitute for filling of intraoral / maxillofacial osseous defects. Indications may include:

- Defects after removal of bone cysts
- Periodontal defects
- Augmentation of alveolar crest (e.g. guided bone regeneration, GBR)
- Defects after root end resection
- Extraction sockets
- Defects after surgical removal of retained teeth
- Sinus floor augmentation
- Defects after removal of autologous bone

Contraindications

calc-i-oss™ should not be used if the patient's health does not allow oral surgery.

calc-i-oss™ should not be applied in case of the presence of any diseases or therapies which are detrimental to the healing of bone (e.g. hyperparathyroidism, osteomalacia, chronic high dose therapy with corticosteroids, liver disease, alcohol use).

Adverse effects

Risks and adverse reactions by the use of calc-i-oss™ correspond to general risks and adverse effects associated with autogenous bone or bone graft substitute grafting procedures, like e.g. post-surgical edema and bleeding, numbness by nerve injury, infection of the soft tissue and/or bone (osteomyelitis), soft tissue ulceration, incomplete bone formation, lack of bone formation, fracture or inflammation of the newly formed bone.

Interactions with other substances

None known.

Warnings

calc-i-oss™ should not be used in inflamed or infected sites.

calc-i-oss™ is sterilized by radiation and is intended for single use. Unused material should be discarded and should not be reused or resterilized.

calc-i-oss™ should not be used if the package is opened or damaged or if the product has expired.

Blood supply to the defect area should not be compromised.

calc-i-oss™ is not intended to bear mechanical loads. The mechanical stability in the region of the defect should be ensured.

Dosage form

calc-i-oss™ is available in the granule sizes 315 – 500 µm, 500 – 1000 µm and 1000 – 1600 µm.

Application

calc-i-oss™ is restricted to sale to or on the order of licensed dentists.

The user should be trained in surgical techniques for dental tissue regeneration and should be experienced in the use of bone graft substitutes.

calc-i-oss™ is delivered sterile and should be opened in an aseptic environment. The material should be used immediately after opening.

The implantation site should be free of infections and soft or granulation tissue. An adequate amount of viable bone should be present at the implantation site.

calc-i-oss™ should be mixed with the patient's blood or blood preparations (e.g. platelet-rich plasma), sterile saline or autogenous bone prior to application.

The underlying bone should be refreshed to generate bleeding. Complete defect filling with bone graft substitute should be attempted. The bone graft substitute should be in intimate contact with the bony defect walls. Do not overfill the defect.

calc-i-oss™ may be used with dental membranes, e.g. to contain the material in the defect, for guided bone regeneration or guided tissue regeneration.

Primary closure should be attempted. The soft tissue coverage should be complete and tension-free.

Smoking diminishes the capacity for bone healing. Inform the patients that they are highly advised to abstain from smoking until completion of bone healing.

Properties

calc-i-oss™ consists of phase-pure β -tricalcium phosphate. The material is osteoconductive and is resorbed through physiological solution in body fluid and direct cellular attack. The resorption will take place mostly parallel to bone regeneration. Depending on the regeneration potential of the tissue, calc-i-oss™ will be completely resorbed within 5 to 15 months.

calc-i-oss™ is solely a bone graft substitute. The application of the material does not guarantee a treatment success. The success of the therapy depends on many factors such as surgical technique and habits, age and the regeneration potential of the patient.

Shelf life

The product should not be used after the imprinted expiry date.

Storage instructions

Store the product in the original packaging under dry conditions.

Release date of the instructions for use

2011-08-08 (original German version)

Manufacturer




Degradable Solutions AG
Wagistrasse 23
8952 Schlieren / Zurich / Switzerland

A company of:

Sunstar Suisse S.A.
Route de Pallatex 15
1163 Etoy / Switzerland

Distribution

Degradable Solutions AG
Wagistrasse 23
8952 Schlieren / Zurich / Switzerland

	Notified body
	For single use only
	Batch number
	Expiration date
	Refer to instruction leaflet
	Sterilisation by gamma irradiation

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