### Description / composition

calc-i-oss™CRYSTAL is a bone graft substitute made from biphasic calcium phosphate (a compound of 60% hydroxyapatite and 40% β-tricalcium phosphate). It is available as granules (granule size 450 – 1000 µm).

calc-i-oss™CRYSTAL is biocompatible and osteoconductive. In the defect, the granules form a macroporous structure that fosters formation of autologous tissue. The granule's high microporosity allows for taking up considerable amounts of blood. With the uptake of blood, factors that are beneficial for healing are made available.

calc-i-oss™CRYSTAL is partially resorbable and remains integrated into bone, preserving the volume.

calc-i-oss™CRYSTAL is synthetic and does not contain animal or human-derived substances.

calc-i-oss™CRYSTAL is sintered at over 1000 °C. The granulation and sinter process produces round, porous granules of high stability with only marginal inclination towards fine particle attrition.

### Indications

calc-i-oss™CRYSTAL is a bone graft substitute for filling of intraoral / maxillofacial osseous defects. Indications may include:

- Defects after removal of bone cysts
- Periodontal defects
- Augmentation of alveolar crest (e.g. guided bone regeneration, GBR)
- Defects after root end resection
- Extraction sockets
- Defects after surgical removal of retained teeth
- Sinus floor augmentation
- Defects after removal of autologous bone

### Contraindications

calc-i-oss™CRYSTAL should not be used if the patient's health does not allow oral surgery.

calc-i-oss™CRYSTAL should not be applied in case of the presence of any diseases or therapies which are detrimental to the healing of bone (e.g. hyperparathyroidism, osteomalacia, chronic high dose therapy with corticosteroids, liver disease, alcohol use).

### Adverse effects

Risks and adverse reactions by the use of calc-i-oss™CRYSTAL correspond to general risks and adverse effects associated with autogenous bone or bone graft substitute grafting procedures, like e.g. post-surgical edema and bleeding, numbness by nerve injury, infection of the soft tissue and/or bone (osteomyelitis), soft tissue ulceration, incomplete bone formation, lack of bone formation, fracture or inflammation of the newly formed bone.

### Interactions with other substances

None known.

### Warnings

calc-i-oss™CRYSTAL should not be used in inflamed or infected sites.

calc-i-oss™CRYSTAL is sterilized by radiation and is intended for single use. Unused material should be discarded and should not be reused or resterilized.

calc-i-oss™CRYSTAL should not be used if the package is opened or damaged or if the product has expired.

Blood supply to the defect area should not be compromised.

calc-i-oss™CRYSTAL is not intended to bear mechanical loads. The mechanical stability in the region of the defect should be ensured.

Since the granules are only partially resorbable, calc-i-oss™CRYSTAL should not be implanted into the surrounding soft tissue.

### Application

calc-i-oss™CRYSTAL is restricted to sale to or on the order of licensed dentists.

The user should be trained in surgical techniques for dental tissue regeneration and should be experienced in the use of bone graft substitutes.

calc-i-oss™CRYSTAL is delivered sterile and should be opened in an aseptic environment. The material should be used immediately after opening.

The implantation site should be free of infection and soft or granulation tissue. An adequate amount of viable bone should be present at the implantation site.

calc-i-oss™CRYSTAL should be mixed with the patient's blood or blood preparations (e.g. platelet-rich plasma), sterile saline or autogenous bone prior to application.

The underlying bone should be refreshed to generate bleeding. Complete defect filling with bone graft substitute should be attempted. The bone graft substitute should be in intimate contact with the bony defect walls. Do not overfill the defect.

calc-i-oss™CRYSTAL may be used with dental membranes, e.g. to contain the material in the defect, for guided bone regeneration or guided tissue regeneration.

Primary closure should be attempted. The soft tissue coverage should be complete and tension-free.

Smoking diminishes the capacity for bone healing. Inform the patients that they are highly advised to abstain from smoking until completion of bone healing.

### Properties

calc-i-oss™CRYSTAL consists of biphasic calcium phosphate. Biphasic calcium phosphate is partially resorbable by cellular mechanisms and direct dissolution. A part of the biphasic calcium phosphate remains embedded in autogenous tissue for prolonged periods of time in order to preserve the volume.

calc-i-oss™CRYSTAL is solely a bone graft substitute. The application of the material does not guarantee a treatment success. The success of the therapy depends on many factors such as surgical technique and habits, age and the regeneration potential of the patient.

### Shelf life

The product should not be used after the imprinted expiry date.

### Storage instructions

Store the product in the original packaging under dry conditions.

### Release date of the instructions for use

2011-08-16 (original German version)

### Manufacturer


Degradable Solutions AG  
Wagistrasse 23  
8952 Schlieren / Zurich / Switzerland

A company of:

Sunstar Suisse S.A.  
Route de Pallatex 15  
1163 Etoy / Switzerland

### Distribution

Degradable Solutions AG  
Wagistrasse 23  
8952 Schlieren / Zurich / Switzerland

	Notified body
	For single use only
	Batch number
	Expiration date
	Refer to instruction leaflet
	Sterilisation by gamma irradiation

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