SUNSTAR
GUIDOR[®]
 bioresorbable
matrix barrier

English

Bioresorbable matrix barrier
 Aid in bone regeneration and augmentation
For Professional Use Only - For Dental Use Only

INDICATIONS FOR USE

The use of the GUIDOR[®] bioresorbable matrix barrier (GUIDOR matrix barrier) to aid in bone regeneration and augmentation in the oral cavity should be limited to defects and concavities within skeletal contours and to defects or situations where moderate increase of bone volume beyond the skeletal contours is desirable. In all cases, appropriate space making support should be used.

CONTRAINDICATIONS

The GUIDOR matrix barrier is contraindicated in those situations where general oral surgery should not be performed. Currently there are no known additional contraindications to the use of the GUIDOR matrix barrier.

PRECAUTIONS

- The GUIDOR matrix barrier is not intended for use in defects other than those stated under INDICATIONS FOR USE.
- The GUIDOR matrix barrier has not been clinically tested:
 - In patients with extra large defects;
 - For extensive bone augmentation;
 - For use in the treatment of failing implants;
 - In immunocompromised patients (diabetes, chemotherapy, irradiation, infection with HIV);
 - In children, pregnant or nursing women;
 - The GUIDOR matrix barrier cannot be reused or resterilized. Resterilization will damage the GUIDOR matrix barrier.

ADVERSE REACTIONS


Possible complications following any oral surgery include thermal sensitivity, flap sloughing, some loss of crest bone height, abscess formation, infection, pain, and complications associated with the use of anesthesia. As with any type of surgical therapy, the patient may experience minor discomfort for a few days. Any serious incident that has occurred in relation to the GUIDOR matrix barrier should be reported to the European authorized representative and the competent authority of the country where the dental professional and/or patient is established.

MATERIAL DESCRIPTION

The GUIDOR matrix barrier is made of an off-white, bioresorbable material consisting of poly-D, L-lactide, poly-L-lactide, and acetyl-tri-butylcitrate (ATBC) a citric acid ester.

At storage temperatures 2°C to 8°C, the GUIDOR matrix barrier is stiff and brittle and may fracture if manipulated. It reaches room temperature after 15 minutes in the foil pouch. At body temperature, the material becomes malleable within seconds.

DESIGN DESCRIPTION

The GUIDOR matrix barrier consists of:

- An external layer (A) facing the mucosal tissue and
- An internal layer (B) facing the jaw bone.
- The two layers are separated by many spaces (C) to form an interspace into which tissue can migrate. The external layer has rectangular perforations (D) (400 – 500 per cm²). The internal layer has circular perforations (E) (4000 – 5000 per cm²). The internal layer includes outer spaces (F) on the jaw bone side of the GUIDOR matrix barrier.

SPACE MAKING

Appropriate space making support for the GUIDOR matrix barrier must be used to allow for bone regeneration. The GUIDOR matrix barrier material is malleable to facilitate application and adaptation of the GUIDOR matrix barrier toward the adjoining bone surface. Without support i.e., sufficient amount of grafting materials), the GUIDOR matrix barrier material may collapse into the defect. Their use for defects with highly favorable morphology i.e., very narrow defects within the skeletal contour, is left to the judgment of the practitioner.

INSTRUCTIONS FOR USE

PREPARATION OF SITE

- Use a flap management technique that ensures total defect coverage.
- Place the primary incision and releasing incisions over intact bone next to the area intended for treatment.
- Remove all granulation tissue from the bone surface.
- In case of implant surgery, follow the manufacturer's instructions for preparation of the implant site and implant insertion.
- Make penetrations through the cortical bone walls into the bone marrow to provide optimal access to bone forming cells and to promote blood clot formation, at the discretion of the practitioner.

PLACEMENT OF THE GUIDOR matrix barrier

- Maintain a sterile field throughout the procedure.
- Make sure the GUIDOR matrix barrier has reached room temperature before manipulating or trimming. Do not attempt to warm the GUIDOR matrix barrier in your hands.
- After the GUIDOR matrix barrier reaches room temperature, remove it from the foil pouch and open plastic tray.
- Identify the external layer of the GUIDOR matrix barrier. In the plastic tray the external layer faces the lid. The external layer can be identified by the "S" marked on this side and is also marked by multiple 2 mm square notches in the margin of this layer.
- Remove the GUIDOR matrix barrier from the plastic tray by lifting the GUIDOR matrix barrier with a flat tweezer.
- To trim, hold the GUIDOR matrix barrier gently between the fingers and cut with sharp scissors.
- Place grafting material for appropriate space making support. In case of an exposed implant surface, place a minimum 2 mm thick layer of grafting material.
- The external layer with the "S" mark should face the mucosal tissue after placement of the GUIDOR matrix barrier.

- When shaping the GUIDOR matrix barrier, ensure the size is sufficient to cover the defect and extend at least 3 mm onto intact bone to prevent the GUIDOR matrix barrier from collapsing into the defect.
- Avoid applying excessive pressure on the GUIDOR matrix barrier throughout the procedure to avoid damage to the GUIDOR matrix barrier.

WOULD CLOSURE

- Make every attempt to achieve primary closure to promote primary healing.
- Avoid tension in the repositioned flap.
- In cases where the GUIDOR matrix barrier is left exposed, it is recommended not to use wound dressings, or liquid barriers, that may adhere to the membrane or obstruct the rectangular holes of the external layer. Movement of the wound dressing could dislodge or damage the membrane and impede healing.

GUIDOR matrix barrier EXPOSURE

Exposed GUIDOR matrix barrier material does not need to be trimmed. It will disappear 6 to 8 weeks after surgery due to combined resorption and attrition.

USE OF ANTIOTIBIICS

Antibiotic therapy may be provided at the discretion of the practitioner and should follow current standards of care. In cases of infection or abscess, removal of the GUIDOR matrix barrier may be indicated depending on the severity of the complication.

POST-SURGICAL CARE

Special instructions are needed for the surgically treated area(s). Avoid pressure by palpation or a prosthesis on the treated area during the healing period. Following surgery, the patient should refrain from mechanical cleaning techniques that are close to the treated area for a minimum of 4 weeks after surgery. During this period mouth rinsing with an antimicrobial agent such as chlorhexidine is recommended.

RE-ENTRY

The area(s) treated for bone regeneration or augmentation should not be uncovered for approximately 4 to 6 months after placement of the GUIDOR matrix barrier. Abutment connection surgery to previously placed implants may still be performed at an earlier date if disturbance to the area treated for bone regeneration or augmentation can be avoided.

Remnants of the GUIDOR matrix barrier material may be present at the time of re-entry, usually within the soft tissue flap. These remnants have a typically liquefied and/or granular appearance. No special measures need to be undertaken, since the remnants will resorb and disappear.

PRECLINICAL STUDIES

The theoretical basis for the function of the GUIDOR matrix barrier has been evaluated in vivo and histologically evaluated following treatment of recession type defects and intrabony defects in monkeys. The histological analysis showed ingrowth of gingival connective tissue through the proportionally large perforations of the external layer of the GUIDOR matrix barrier to fill out the interspace, thereby preventing or minimizing epithelial downgrowth. The relatively small perforations of the internal layer resulted in an adapted delay of ingrowth of connective layer. Integration of the GUIDOR matrix barrier with the soft tissue flap during the initial healing also minimized gingival recession and GUIDOR matrix exposure. The results showed that, when implanted and integrated with the tissues, the design of the GUIDOR matrix barrier was fully preserved for at least 6 weeks after surgery. By 3 months the resorption process was initiated and the GUIDOR matrix barrier had begun to fragment. Bioreosorption took up to 12 months to be completed. New attachment and new supporting bone had formed.

The function of the GUIDOR matrix barrier in bone regeneration and augmentation has been studied clinically and histologically in rats and rabbits. Healing and integration of the GUIDOR matrix barrier with the surrounding tissues were unevenful. Re-entry surgery of bone augmentation in rabbits showed that significantly more new bone had formed when the GUIDOR matrix barrier was used in combination with supporting autologous bone, compared to when autologous bone was used alone.

PACKAGING

Contains one GUIDOR matrix barrier. The GUIDOR matrix barrier comes sterilized (E-Beam) in a transparent plastic tray designed to protect the GUIDOR matrix barrier from physical damage. The plastic tray is packed in a foil pouch, which maintains the sterility of the GUIDOR matrix barrier. Do not use if package is damaged. Store between 2°C to 8°C.

The GUIDOR name and logo are registered trademarks of Sunstar Suisse SA in Europe and in other countries.

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